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Instructions for Use TPP Tissue Culture Flasks

TPP tissue culture flasks are intended for cell/tissue culture in manual use. The growth surface - but not the side wall of the flask – is opto-mechanically activated for optimal cell adhesion and growth. The angled neck reduces the risk of medium contacting the inside of the screw cap by spillage and allows excellent access for cell scrapers and serological pipettes.

TPP tissue culture flasks are available with filter or VENT screw caps.

The tissue culture flask is intended for single use only.

Note:

Follow national regulations when handling biological material, use appropriate protective clothing

Take note of the regulations of aseptic working during the work process.

Instructions:

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Check the expiry date (EXP) on label and packaging. Only use products with valid EXP.*
- *Open the packaging in a sterile environment and remove a complete system for use.*
- *Open the flask and fill it with the medium and inoculum according to your laboratory routine.*

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Filter Screw Cap:

- For the continuous gas exchange, close the filter screw cap with integrated hydrophobic PTFE filter membrane 0.22 µm tightly (picture 1). If the membrane gets wet, the gas exchange reduces for a short period.

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
VENT Screw Cap:

- For a gas exchange, turn the VENT cap into the ventilation position at 12 o'clock. You will hear/notice a click.
 - Visual check:
Rectangle at 12 o'clock (position 1)
- During incubation of cells, the rectangle must be in position 1 / 12 o'clock.
- To interrupt the gas exchange, turn the cap a 1/4 turn clockwise over the aeration position. The flask is now tightly close there is no gas exchange.
 - Visual check:
Rectangle at 3 o'clock. (position 2)
- For a transfer of the flask within the laboratory (laminar flow/incubator), TPP recommends to close the VENT screw cap tightly.

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
Advice and Recommendations:

- *For labelling there are bilateral inscription fields*
- *The graduation serves as a reference for the filling quantity and is not an absolute value.*
- *The stacking rim allows a stable stacking of several flasks on top of each other, also in combination with other TPP flasks of the same size.*
- *For a homogenous cultivation of the cells, levelling of the incubator and continuous supply of medium is essential.*
- *Avoid and/or eliminate vibration sources near the incubator to ensure uniform cell growth.*

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Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch