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Instruction for Use TPP TubeSpin® Bioreactor

The TubeSpin® bioreactor is for high-throughput screening, cultivation and optimization processes of prokaryotic (aerobic bacteria) and eukaryotic (yeasts, algae, HEK, CHO, Sf-9) suspension cells.

The TubeSpin® bioreactor is available with a filter screw cap with 5 or 10 openings over a PTFE filter membrane.

The TubeSpin® bioreactor is for single use only.

Note:

Follow the national regulations for handling biological materials and wear the appropriate protective clothing.

During the work process, be sure to observe the rules of aseptic technique.

Instructions:

- *Inspect packaging and product for integrity. Use only materials in perfect condition.*
- *Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.*
- *Open the package in a sterile environment and remove a complete system for use.*
- *Open the bioreactor and fill it with the medium and inoculum according to your laboratory routine. Please refer to the optimal filling volume (see technical data).*
- *Close the filter screw cap with ventilation holes. Continuous gas exchange occurs through the integrated 0.22 µm hydrophobic membrane.*

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
Bild 1

- *If the gas exchange needs to be reduced, individual openings can be sealed with laboratory grade labels or tape (Figure 1).*


Figure 1

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- **Cultivation:** Shake the cells in appropriate incubation shakers. Observe the cultivation conditions of the cells (temperature, CO₂ and humidity control, shaking frequency).
- **Shakers:** Follow the manufacturer's safety and operating instructions.
- **Centrifugation:** The use of appropriate rotors or centrifuge adapters is recommended. Follow the centrifuge manufacturer's safety instructions.
- The mechanical strength of the tubes during centrifugation is influenced by
 - Shape and material
 - Accuracy of fit in the centrifuge adapter
 - Temperature, centrifugation time, relative centrifugal force (RCF)
 - Chemical and physical properties of the centrifuged material
 - Rotor type: fixed-angle or swing-out rotor
- Test the many influencing factors in advance under routine conditions.
- Store at room temperature, protected from UV light.

Remarks and Recommendations:

- *The volume scale is an indication of the filling volume, not an absolute value.*
- *Alcohol 90% with mechanical stress (rubbing) may dissolve or blur the white color of the print.*

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Storage in Minus Temperatures:

- *The product should only be used for its intended use. TPP does not assume any product warranty in case of damage when used otherwise.*
- *For storage in temperatures below °C transfer the samples into TPP cryo tubes. If centrifuge tube storage in the minus temperature range is desired, TPP recommends testing the product in advance under the selected conditions.*

Technical Data:

Materials

Screw Cap	PE
Membrane	PTFE
Tube	PP

Additional:

Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from www.tpp.ch.

Abmessungen	Measurements	87015	87017	87050	87056
Volumen grad. ml	<i>Volume grad. mL</i>	15	15	50	50
Länge mm	<i>Length mm</i>	120	120	115	115
Durchmesser mm	<i>Diameter mm</i>	17.1	17.1	30	30
Max. RZB x g	<i>Max. RCF x g</i>	15'500	15'500	15'500	15'500
Form	<i>Form</i>	konisch/ <i>conical</i>	rund/ <i>round</i>	konisch/ <i>conical</i>	rund/ <i>round</i>
Optimale Einfüllmenge ml	<i>Optimal filling volume mL</i>	1 – 10	1 – 10	1 – 35	1 – 35
Schüttler: Empf. Orbit / Schütteldurchmesser mm	<i>Shaker: Recom. Orbit / Shaking diameter mm</i>	50	50	50	50
Schüttler: Empf. Geschwindigkeit rpm	<i>Shaker: Recom. speed rpm</i>	250	250	180	180
Schraubkappe "solid"	<i>Screw cap „solid“</i>	--	--	--	--
Schraubkappe "filter"	<i>Screw cap „filter“</i>	--	--	--	--

Abmessungen	Measurements	87450	87600
Volumen grad. ml	<i>Volume grad. mL</i>	450	600
Länge mm	<i>Length mm</i>	145	181
Durchmesser mm	<i>Diameter mm</i>	96	98
Max. RZB x g	<i>Max. RCF x g</i>	3'500	3'500
Form	<i>Form</i>	Konisch / <i>conical</i>	Konisch / <i>conical</i>
Optimale Einfüllmenge ml	<i>Optimal filling volume mL</i>	>300	>400

Schüttler & Schraubkappe	Shaker & Screw Cap	87450	87600
Schüttler: Empf. Orbit / Schütteldurchmesser mm	<i>Shaker: Recom. Orbit / Shaking diameter mm</i>	50	50
Schüttler: Empf. Geschwindigkeit rpm	<i>Shaker: Recom. speed rpm</i>	150	150
Schraubkappe "solid"	<i>Screw cap „solid“</i>	87355	87355
Schraubkappe "filter"	<i>Screw cap „filter“</i>	87356	87356

Öffnungen / Openings Ø mm	87015 / 87017	87050 / 87056	87450 / 87600
A	0.4	0.4	-
B	0.6	0.6	-
C	1.0	1.0	-
D	1.25	1.5	-
E	1.5	2.0	-
10 x each	-	-	4.0