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Instruction for Use TPP Cryo Tube

TPP cryo tubes are made of polypropylene (PP) with external thread and star feet are designed for storage of biological material and samples down to -196 °C. The screw caps with external thread reduce the risk of contamination and cross-contamination.

The cryo tube is for single use.

Note:

Follow the national regulations for handling biological materials and wear the appropriate protective clothing.

During the work process, be sure to observe the rules of aseptic technique.

TPP recommends storing the tubes only in the gas phase of nitrogen (LN₂). If stored in the liquid phase, the LN₂ could penetrate. During the thawing process, rapid evaporation would cause internal pressure to build up. This could result in an explosion of the tube and the release of biologically hazardous substances.

Instructions:

- *Inspect packaging and product for integrity. Use only materials in perfect condition.*
- *Check the expiry date (EXP) on the label and packaging. Use only products with a valid EXP.*
- *Open the package in a sterile environment and remove a complete system for use.*
- *For a firm and secure stand while working, place the cryotubes in the cryotube rack (#99016).*

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Storage in Minus Temperatures:

- *When storing the cryo tubes in the subzero range, first test the many influencing factors under routine conditions.*
- *Temperatures below 0 °C result in reduced elasticity for polypropylene. Do not subject frozen tubes to mechanical stress, such as dropping or hitting it.*
- *Use only "freezing containers" for controlled freezing of cells before transferring the cryo tubes to the gas phase of a nitrogen tank.*

Advice and Recommendations:

- *The volume scale is an indication of the filling volume, not an absolute value.*
- *A labeling area and colored cap inserts allow for sample identification.*
- *An example of a cell freezing protocol can be found on our website at www.tpp.ch.*

Accessories:

- *Cryo Box L, # 99014 / Cryo Box S, # 99015*
- *Cryo Rack # 99016: The star foot in combination with the TPP rack allows one-handed operation.*
- *Colored cap inserts for identification*

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Instructions:

Designed for the storage of cryo tubes and most other similarly designed vials at ultra-low temperatures.

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Additional:

Instructions for use, chemical resistance lists, and quality certificates for each product are available for download from www.tpp.ch.