



Andreas Hettich GmbH & Co KG

Centrifuges and Rotor Compatibility TubeSpin® Bioreactor 450

The following rotors of the centrifuge from Hettich are compatible for the TubeSpin® Bioreactor 450.



Rotanta 460



Rotanta 460R

Hettich® ROTANTA 460/460R Centrifuge

Table top centrifuge 5650/5660



5699-R
Rotor



4449
Adapter



4890
Hanger with cup
cup bio-safe 4883



4895
Hanger without closure
use with adapter



Andreas Hettich GmbH & Co KG

Centrifuges and Rotor Compatibility TubeSpin® Bioreactor 450 and 600

The following rotors of the centrifuge from Hettich are compatible for both TubeSpin® Bioreactor 450 and 600




Roto Silenta 630 RS

Hettich® ROTO SILENTA 630 RS centrifuge

Table top centrifuge 5650/5660



4176
Rotor



SK 11.10-1
Adapter



SK 11.10
Hanger with cup
cup bio-safe 4883



4895
Hanger without closure
use with adapter

For more technical information on centrifuges, rotors and adapters from
 Andreas Hettich GmbH & Co KG
 please visit www.hettichlab.com and/or contact info@hettichlab.com.



Technical Data
TubeSpin® Bioreactor 450 / 600

Measurements	87450	87600
Volume graduation mL	450	600
Length mm	145	183
Diameter mm	96	100
Max. RCF x g	3'500	3'500
Form	conical	conical
Optimal filling volume mL	>300	>400
Shaker: Recom. orbit/shaking throw mm	50	50
Shaker: Recom. speed rpm	150	150
Screw cap „solid“	87355	87355
Screw cap „filter“	87356	87356