



## Eppendorf Centrifuges and Rotor Compatibility TubeSpin® Bioreactor 600

The following rotor of the Eppendorf 5920 R centrifuge is compatible for the TubeSpin® Bioreactor 600.



**5920 R**

Swing out-rotor		S-4x1000	5895 101 003 <sup>[1]</sup>
Bucket	Bucket Lid	Adapter	Capacity
5895 102 000 (set of 4 pcs.)	5820 747 005 (set of 2 pcs.)	5920 701 002 (set of 2 pcs.)	4 x 87600

<sup>[1]</sup> without buckets

For more technical information on centrifuges, rotors and adapters from  
 Eppendorf  
 please visit <https://www.eppendorf.com> and/or contact your local dealer.

### Technical Data TubeSpin® Bioreactor 450 / 600

Measurements	87450	87600
Volume graduation mL	450	600
Length mm	145	183
Diameter mm	96	100
Max. RCF x g	3'500	3'500
Form	conical	conical
Optimal filling volume mL	>300	>400
Shaker: Recom. orbit/shaking throw mm	50	50
Shaker: Recom. speed rpm	150	150
Screw cap „solid“	87355	87355
Screw cap „filter“	87356	87356