




Mounting & Disassembling a TPP Tube Rack



Rack "S" and "L" for centrifuge tubes

- Identification: Alphanumeric coded
- Storage: Space saving, as it can be quickly disassembled
- Holds:
 - Rack "S": 18 x 15 mL and 10 x 50 mL
 - Rack "L": 30 x 15 mL and 20 x 50 mL

| Measurements | 99017 | 99019 |
|---------------------------|------------------------------------|-------|
| Length mm | 209 | 209 |
| Width mm | 98 | 172 |
| Height mm | 62 | 62 |
| Tubes per rack | 18 | 30 |
| Materials | PP | |
| Working Temperature | +121 to -190 °C | |
| Autoclavable | yes | |
| Sterilizing instructions: | Please sterilize disassembled only | |